



# GloMax DISCOVER

A high-performance, easy-to-use multimode reader for luminescence, fluorescence, absorbance, BRET and FRET applications.



# The Perfect Partner for Promega Assays

## Easy-To-Use

Choose from preloaded Promega protocols or customize your own protocol. Export data to a network, cloud, LIMS or any drive desired. Easily perform luminescent, fluorescent, UV-visible absorbance, BRET and FRET measurements.

## Automation-Ready

Integrate into your automated workflow for higher walk-away throughput or into your LIMS data management system.

## Superior Sensitivity

Broader dynamic range, better sensitivity and lower well-to-well cross talk for more usable data.

## Service

Full-service products are available to meet your needs including Installation Qualification (IQ) and Operational Qualification (OQ).

## Minimal Manual Intervention

Automated filter switching allows easy multi-assay detection.

## Cost Effective

Preloaded protocols save time on assay optimization; less hands-on time conserves personnel resources.

## Reporting

Provides the required technical elements of a 21CFR Part 11 compliant system for electronic signature control.


## Intuitive Software Interface

Create custom protocols, easily integrate into automated platforms, and export your data—all with the tap of a screen.

- Choose from preloaded protocols or customize your own
- Set up a plate map
- View data in a heat map display
- Easily export your results



Heat Map Display



Drag and Drop Protocols

# Assays that Integrate with Your Instrument

Choose from preloaded Promega protocols or customize your own using the intuitive drag-and-drop interface. GloMax<sup>®</sup> Discover was designed and developed for Promega assays:

- Making your research experiments easy
- Getting you up and running faster
- Supporting all of your needs in one place (assay, instrument, service)

## Cell Health Assays:

- CellTiter-Glo<sup>®</sup>
- CellTox<sup>™</sup> Green
- Caspase-Glo<sup>®</sup>
- *and many more*

## Cell Signaling & Metabolism Assays:

- ADP-Glo<sup>™</sup>
- Kinase-Glo<sup>®</sup>
- P450-Glo<sup>®</sup>
- *and many more*

## BRET and FRET Assays:

- NanoBRET<sup>™</sup>
- Commercial and Homebrew assays

## Luciferase Reporter Assays:

- Nano-Glo<sup>®</sup>
- ONE-Glo<sup>®</sup>
- Dual-Glo<sup>®</sup> & DLR<sup>™</sup>
- *and many more*



Integrate  
Into Existing  
Workflows

NanoLuc<sup>®</sup>  
Luciferase  
Ready


DLR<sup>ready</sup>

Perform  
Multiplex  
Assays

# Reliable Data Collection

## Luminescence


- Low-noise PMT ensures collected light is not compromised
- 10 to 1,000 times more sensitive than competitor instruments
- No need to dilute samples or manage Gain settings
- 2 to 3 logs more usable data than competing instruments
- Total or filtered luminescence readout for BRET or multi-wavelength distinction
- Lower cross talk for more usable data
- 4 filters included for BRET and filtered luminescence assays such as NanoBRET™, Renilla/YFP BRET, and Chroma-Glo™; Plus 1 customizable position
- Or, customize all 5 luminescence filters to your needs



Tyrosine kinase inhibition was assessed following incubation with Gefitinib, a selective inhibitor of EGFR, using the ADP-Glo™ Kinase Assay.

## Luminescence Module Specifications

Detector	Top-reading, head-on photon counting photomultiplier tube (PMT)
Wavelength Range	350–650nm
Detection Limit	$3 \times 10^{-21}$ moles of luciferase
Linear Dynamic Range	9 logs
Cross talk	Less than $3 \times 10^{-5}$ (white, 96-well plate, Corning 3912)
Filtered Luminescence	Built-in filter paddle



Multiplexed cytotoxicity and viability assays were performed on digitonin-treated K562 cells.

Study cell viability, apoptosis and cytotoxicity

## Fluorescence

- Uses powerful wavelength-matched light-emitting diodes (LED) as excitation sources
- LEDs closely match excitation profiles of commonly used fluorescent molecules to ensure high sensitivity
- Reduces nonspecific light leakage, a common problem when using broad-spectrum light sources
- 5 excitation and emission filters included for common fluorescence assays; Plus 1 customizable excitation and emission position
- Or, customize two fluorescence excitation positions or all six emission positions to your needs

### Fluorescence Module Specifications

Detector	PIN-photodiode
Light Source	Wavelength-matched LED
Read Position	Top reading
Wavelengths selection	Filter module with 5 standard excitation and emission filters. Empty filter positions are available for custom excitation and emission wavelengths.
Wavelengths included	UV (Ex: 365nm, Em: 415–445nm) Blue (Ex: 475nm, Em: 500–550nm) Green (Ex: 525nm, Em: 580–640nm) Red (Ex: 625nm, Em: 660–720nm) AFC (Ex: 405nm, Em: 495–505nm)
Detection Limit	2fmol fluorescein/200µl
Linear Dynamic Range	>6 logs (assay dependent)

## UV-Visible Absorbance

- Comes factory-installed with 9 UV-Visible filters for common colorimetric assays
- Ideal for ELISA, protein, and nucleic acid assays

### UV-Visible Absorbance Module Specifications

Detector	Top-reading, head-on photon counting photomultiplier tube (PMT)
Light Source	Xenon flash lamp
Spectra Range	200–600nm
Filter Wheel	9 absorbance filters included
Wavelengths included	230, 260, 280, 320, 405, 450, 490, 560 and 600nm 10nm bandwidth
Detection Limit	0.1 O.D.
Dynamic Range	0–4.0 O.D.



Integrate Into Your Automated Workflow

*GloMax® Discover is Compatible with  
Third Party Automation Control*

## GloMax<sup>®</sup> Discover System Specifications

### Instrument Specifications

Detection Modes	Luminescence, Fluorescence Intensity, UV-Visible Absorbance, Filtered Luminescence, BRET and FRET
Sample Format	6, 12, 24, 48, 96 and 384-well plates
Tablet PC	Windows <sup>®</sup> 8 Tablet PC with USB connection to instrument
Read Speed	96-well plate: less than 1 minute; 384-well plate: less than 3 minutes
Automation Compatibility	Hardware and software accessible to 3rd party hardware/software control, including LIMS.
Dimensions	18.5in (46.9cm) width, 19.5in (46.5cm) depth, 10in (25.4cm) height without Tablet and PC mounted, 16in (35.5cm) height with Tablet and PC mounted
Weight	Approximately 60lbs. (27kg)
Power Requirements	120VA, 50/60 HZ
Warranty	1-year parts and labor warranty included
Regulatory Compliance	For Research Use Only. Not for use in Diagnostic Procedures

### Injector System Specifications

Number of injectors (optional)	Two injectors
Dispense Volume Range	5–200µl in 1µl increments
Plate Compatibility	6, 12, 24, 48, and 96-well plates
Injection Speed	20–500µl per second
Waste Collection Tray Volume	Approximately 50ml
Void Volume	500µl

### Shaker System Specifications

Shaking Pattern	Linear or orbital patterns
Shaking Intensity	100–500 cycles/minute

### Heating System Specifications

Temperature Range	Ambient +5°C to 45°C
-------------------	----------------------

### PC Requirements

Operating System	Windows <sup>®</sup> 8, 64 bit, with USB connection to instrument
RAM	2GB
Disk Space	64GB
Other Requirements	Microsoft.NET 4.0 Framework (full version), Microsoft Excel 2013 Professional



## Ordering Information

GloMax® Discover System <i>(comes with U.S. Instrument power cord and U.S. Tablet PC power cord)</i>	Cat.# GM3000
GloMax® Discover Dual Injectors with Pump Station	Cat.# GM3030
GloMax® Discover Luminescence Filter Paddle	Cat.# GM3011
GloMax® Discover Fluorescence Filter Paddle	Cat.# GM3012
GloMax® Discover Inner and Outer Injector Tubing, Natural	Cat.# GM3014
GloMax® Discover Inner and Outer Injector Tubing, Black	Cat.# GM3015
Installation Qualification (IQ)	Cat.# SA1104
Operational Qualification (OQ)	Cat.# SA1105
Installation and Operational Qualification (IQ/OQ)	Cat.# SA1106
Standard Service Agreement	Cat.# SA4000
Basic Service Agreement	Cat.# SA4010
Preventive Maintenance	Cat.# SA4030

Products may be covered by pending or issued patents or may have certain limitations. Please visit our Web site for more information.

Caspase-Glo, CellTiter-Glo, Dual-Glo, GloMax, Kinase-Glo, Nano-Glo, ONE-Glo and P450-Glo are registered trademarks of Promega Corporation. ADP-Glo, CellTox, Chroma-Glo, DLR and NanoBRET are trademarks of Promega Corporation.

Windows is a registered trademark of Microsoft Corporation.

