



eisco

INDUSTRIAL

CERAMICS

- Evaporating Dishes
- Crucibles
- Funnels
- Mortar & Pestles



Evaporating Dishes - Porcelain, Flat Form ¹

- Porcelain, flat form with spout
- Glazed inside and outside except for rim

Code	Capacity	O.D.	Height	Pack Size	Price/Pack ₹
CH0070A	35 ml	65 mm	16 mm	10	831.00
CH0070B	50 ml	75 mm	20 mm	10	929.00
CH0070C	100 ml	100 mm	23 mm	10	1552.00
CH0070D	200 ml	120 mm	34 mm	10	1858.00
CH0070E	300 ml	150 mm	29 mm	10	2799.00



Evaporating Dishes - Porcelain, Round Form ²

- Porcelain, round form with spout
- Fully glazed except outside base

Code	Capacity	O.D.	Height	Pack Size	Price/Pack ₹
CH0071A	35 ml	60 mm	24 mm	10	953.00
CH0071B	70 ml	75 mm	30 mm	10	1137.00
CH0071C	80 ml	80 mm	30 mm	10	1100.00
CH0071D	120 ml	94 mm	42 mm	10	1222.00
CH0071E	150 ml	100 mm	42 mm	10	1430.00
CH0071F	250 ml	115 mm	45 mm	10	2799.00
CH0071G	385 ml	145 mm	48 mm	10	3324.00
CH0071H	525 ml	162 mm	51 mm	1	560.00



Evaporating Dishes - Silica Round Bottom ³

- Round bottom, with spout

Code	Capacity	Dia.	Height	Price/Piece ₹
CH0078A	20 ml	55 mm	23 mm	1177.00
CH0078B	45 ml	75 mm	27 mm	1700.00
CH0078C	70 ml	87 mm	33 mm	2092.00
CH0078D	85 ml	94 mm	36 mm	2420.00
CH0078E	100 ml	100 mm	45 mm	2746.00
CH0078F	200 ml	113 mm	47 mm	4381.00





Crucibles - Porcelain, Tall, With Lid 1

- Porcelain, tall form with lid
- Glazed inside and outside except base



Code	Capacity	Height	Pack Size	Price/Pack ₹
CH0074A	10 ml	26 mm	10	746.00
CH0074B	15 ml	29 mm	10	929.00
CH0074D	30 ml	37 mm	10	1014.00
CH0074E	40 ml	40 mm	10	1100.00
CH0074F	50 ml	43 mm	10	1577.00
CH0074G	100 ml	54 mm	10	2078.00
CH0074H	250 ml	72 mm	1	299.00

Crucibles - Porcelain, Squat, With Lid 2

- Porcelain, squat form with lid
- Glazed inside and outside except base



Code	Capacity	Height	Pack Size	Price/Pack ₹
CH0076A	8 ml	20 mm	10	648.00
CH0076B	17 ml	26 mm	10	929.00
CH0076C	30 ml	30 mm	10	1112.00
CH0076D	50 ml	37 mm	10	1491.00
CH0076E	100 ml	46 mm	10	2237.00
CH0076F	150 ml	52 mm	10	2787.00
CH0076G	250 ml	60 mm	1	458.00

Mortar and Pestle

- Full porcelain construction
- Heavy duty construction resistant to chips and cracks
- Unglazed grinding surface for precise crushing and grinding
- Oversized pestles for efficient grinding



Code	Outer Dia.	Capacity	Price/Piece ₹
CH0540A 3	80 mm	60 ml	156.00
CH0540B	105 mm	150 ml	327.00
CH0540C	125 mm	275 ml	359.00
CH0540D	155 mm	550 ml	653.00
CH0540E	180 mm	900 ml	882.00
CH0540F	200 mm	1500 ml	1510.00
CH0540G	230 mm	2000 ml	1764.00



Buchner Funnels - Porcelain 1

- Buchner funnel with straight sides above the perforated plate
- Glazed inside and outside except the rim
- Thick stem with angled cut

Code	Dia.	Recommended filler Paper Size	Price/Piece ₹
CH0462A	3.5 cm	27 mm	183.00
CH0462B	5 cm	40 mm	229.00
CH0462C	6 cm	50 mm	241.00
CH0462D	7.5 cm	70 mm	339.00
CH0462E	10 cm	90 mm	491.00
CH0462F	12.5 cm	110 mm	692.00
CH0462G	15 cm	130 mm	1196.00
CH0462H	20 cm	175 mm	3354.00
CH0462I	25 cm	230 mm	4838.00
CH0462J	30 cm	260 mm	9154.00



Check out
our website
www.eiscoindustrial.in