

# eisco

INDUSTRIAL

## LAB SAFETY

- Gloves
- Safety Goggles
- Lab Coats
- Aprons





### Disposable Gloves - Nitrile

- Latex free
- Textured finish on fingers for better grip
- Ideal choice for chemical lab and use during dissection
- Non-sterile, disposable
- Supplied in Pack of 100

Code	Specs	Price/Pack ₹
CH0612H	Small	2689.00
CH0612I	Medium	2689.00
CH0612J <sup>1</sup>	Large	2689.00



### Hand Protectors

- Heat & cold resistant silicon rubber hand protector
- Allows easy handling of heated vessels like flask, bottles & beakers etc
- Studded surface is used for firm grip
- Pockets for fingers & thumb, ambidextrous
- Tested for use up to 250°C


Code	Specs	Pack Size	Price/Pack ₹
CH0616HP <sup>2</sup>	Hand Protectors	10	18578.00



### Safety Glasses - Clear

- Designed to fit over most prescription glasses
- Side vents prevent glasses from fogging
- Lens coating provides UV protection
- Light weight polycarbonate, impact resistant lenses
- CE marked, conforms to ANSI Z87.1, EN166F side shields protection

Code	Specs	Pack Size	Price/Pack ₹
CH0619A <sup>3</sup>	Safety Goggles - Clear	10	3911.00



### Safety Goggles - Basic

- Designed to fit over most prescription glasses
- Soft, flexible vinyl frame provides a comfortable and secure fit
- Freely adjustable, soft elastic head strap
- Direct side vents prevents fogging
- Designed to protect eyes from chemical splashes and high speed impact (120m/s)
- CE marked, conforms to ANSI Z87.1, EN166 B

Code	Specs	Price/Piece ₹
G103 <sup>4</sup>	Safety Goggles - Basic	452.00





### Lab Coats - White Polyester 1

- Modern styling, white polyester/cotton drill
- Concealed press stud fastenings, backvent
- Three large pockets for plenty of storage

Code	Size	Pack Size	Price/Pack ₹
CH0622C	Small	10	<b>23222.00</b>
CH0622A	Medium	10	<b>23222.00</b>
CH0622B	Large	10	<b>23222.00</b>
CH0622D	Extra Large	10	<b>23222.00</b>



### Vinyl Aprons 2

- Transparent vinyl aprons are suitable for general lab use
- Protects clothes from acid, oils and grease
- Halter design for unrestricted arm movement
- Reinforced at strain points to resist stretching

Code	Specs	Pack Size	Price/Pack ₹
VNAPCH	Small	50	<b>17111.00</b>
VNAPAD	Large	50	<b>21389.00</b>



### Rubber Lab Aprons - Black

- Suitable for general lab work or for protection during cleaning, repairing or painting
- Waterproof and resistant of chemicals, oils and grease
- Halter design for unrestricted arm movement
- Flexible and light weight

Code	Specs	Pack Size	Price/Pack ₹
RBAPCH	Small	25	<b>13200.00</b>
RBAPMD 3	Medium	25	<b>13200.00</b>
RBAPAD	Large	25	<b>16500.00</b>

### Size Chart

Code	Length (cm)	Chest (cm)	Arm Length (cm)
CH0622A	100	60	46
CH0622B	105	64	50
CH0622C	95	58	42
CH0622D	110	67	54

Code	Width (in)	Overall Length (in)
RBAPCH	24"	30"
RBAPMD	27"	36"
RBAPAD	27"	42"
VNAPCH	27"	36"
VNAPAD	36"	45"

