



Neoprene Tubing - For Burners

- Neoprene with reinforced ends, 5 mm inner dia.
- Fits burner inlet up to 8 mm dia.

Code	Specs	Pack Size	Price/Pack ₹
CH0100D ①	45 cm	10	2664.00
CH0100E	60 cm	10	3269.00
CH0100F	90 cm	10	3948.00



Manual Pipette Fillers

- Manual bulb filler for pipetting, suitable for noxious solutions
- For use with microbial suspensions and other harmful or unpleasant solutions
- Comprised of a bulb and three pinch valves controlling influx and efflux of solutions
- Bulb dia. 54 mm

Code	Specs	Pack Size	Price/Pack ₹
CH0580A ②	Pipette Filler - Rubber	10	2689.00
CH0580D ③	Pipette Filler - Silicon	10	6661.00



Pipette Bulbs - Rubber

- Pear shaped, for use with pipettes

Code	Capacity	Pack Size	Price/Pack ₹
CH0230A ④	5 ml	10	403.00
CH0230B	10 ml	10	440.00
CH0230C	15 ml	10	501.00
CH0230D ⑤	25 ml	10	917.00
CH0230E	50 ml	10	1582.00
CH0230F	20 ml	10	660.00
CH0230G	100 ml	10	2628.00



Pipette Bulbs - Silicon

- Pear shaped, for use with pipettes

Code	Capacity	Pack Size	Price/Pack ₹
CH0230A-SI	5 ml	10	721.00
CH0230B-SI	10 ml	10	990.00
CH0230C-SI	15 ml	10	1296.00
CH0230D-SI	25 ml	10	2713.00
CH0230E-SI ¹	50 ml	10	4913.00
CH0230F-SI	20 ml	10	1381.00
CH0230G-SI	100 ml	10	6588.00



Dropper Bulbs

- For use with dropper pipettes
- Capacity 1 ml approx

Code	Specs	Pack Size	Price/Pack ₹
CH0232A	Bulbs - Rubber	10	177.00
CH0232AA ²	Bulbs - Silicon	10	367.00



Dropper Bulbs - Heavy Duty

- Heavy weight, for use with pipettes & medicine droppers
- Dia. 5 - 6 mm
- Capacity 2 ml approx

Code	Specs	Pack Size	Price/Pack ₹
CH0232B	Bulbs - Rubber	10	367.00
CH0232BB ³	Bulbs - Silicon	10	746.00



Come chat with us
Live on www.eiscoindustrial.in



Rubber Stoppers - Solid 1

- Solid, made from superior quality rubber
- Resistant to acids, alkalis and ammonia

Code	Bottom	Top	Length	Pack	Price/Pack ₹
CH0314A01	6 mm	9 mm	16 mm	10	54.00
CH0314A02	8 mm	10.5 mm	20 mm	10	61.00
CH0314A	9 mm	11.5 mm	20 mm	10	67.00
CH0314B	11 mm	14 mm	24 mm	10	80.00
CH0314C	13 mm	16 mm	24 mm	10	91.00
CH0314D	15 mm	18 mm	24 mm	10	116.00
CH0314E	17 mm	20 mm	26 mm	10	128.00
CH0314F	18 mm	21 mm	26 mm	10	152.00
CH0314G	19 mm	22 mm	28 mm	10	189.00
CH0310H	21 mm	24 mm	28 mm	10	250.00
CH0314I	23 mm	26 mm	28 mm	10	293.00
CH0314J	25 mm	28 mm	28 mm	10	336.00
CH0314K	27 mm	31 mm	32 mm	10	372.00
CH0314L	29 mm	33 mm	32 mm	10	416.00
CH0314M	31 mm	36 mm	35 mm	10	599.00
CH0314N	33 mm	38 mm	38 mm	10	739.00
CH0314O	35 mm	45 mm	36 mm	10	910.00
CH0314P	38 mm	42 mm	40 mm	10	1044.00
CH0314Q	40 mm	49 mm	40 mm	10	1374.00
CH0314R	44 mm	58 mm	45 mm	10	1704.00
CH0314S	48 mm	54 mm	33 mm	10	2640.00
CH0314T	51 mm	60 mm	57 mm	10	2878.00
CH0314X	75 mm	90 mm	39 mm	10	4987.00



1

Rubber Stoppers - Solid, pack of 50

- Solid, made from superior rubber
- Resistant to acids, alkalis and ammonia
- 7 pieces each 8, 9, 11, 13, 15 mm bottom dia.
- 3 pieces each 17, 21, 23, 25, 29 mm bottom dia.
- Pack of 50

Code	Specs	Pack Size	Price/Pack ₹
CH0315 2	Rubber Stoppers - Solid	50	948.00



2



Rubber Stoppers - Solid, Pack of 58 ¹

- Solid, made from superior rubber
- Resistant to acids, alkalis and ammonia
- 8 pieces each 8, 9, 11, 13, 15, 17 mm bottom dia.
- 3 pieces each 21, 23 mm bottom dia.
- 2 pieces each 25, 29 mm bottom dia.
- Pack of 58

Code	Specs	Pack Size	Price/Pack ₹
CH0318	Rubber Stoppers Solid	58	948.00



Rubber Stoppers - One Hole and Two Hole

- Made from superior rubber
- Resistant to acids, alkalis and ammonia
- Available with precision bore 1 hole or 2 holes

1 Hole	2 Hole	Bottom	Top	Length	Pack	Price/Pack ₹	Price/Pack ₹
②	③					1 Hole	2 Hole
CH0314A011H	-	6 mm	9 mm	16 mm	10	78.00	-
CH0314A021H	-	8 mm	10.5 mm	20 mm	10	78.00	-
CH0314A1H	CH0314A2H	9 mm	11.5 mm	20 mm	10	92.00	112.00
CH0314B1H	CH0314B2H	11 mm	14 mm	24 mm	10	107.00	127.00
CH0314C1H	CH0314C2H	13 mm	16 mm	24 mm	10	128.00	141.00
CH0314D1H	CH0314D2H	15 mm	18 mm	24 mm	10	171.00	196.00
CH0314E1H	CH0314E2H	17 mm	20 mm	26 mm	10	171.00	196.00
CH0314F1H	CH0314F2H	18 mm	21 mm	26 mm	10	196.00	213.00
CH0314G1H	CH0314G2H	19 mm	22 mm	28 mm	10	232.00	262.00
CH0314H1H	CH0314H2H	21 mm	24 mm	28 mm	10	293.00	336.00
CH0314I1H	CH0314I2H	23 mm	26 mm	28 mm	10	354.00	403.00
CH0314J1H	CH0314J2H	25 mm	28 mm	28 mm	10	384.00	446.00
CH0314K1H	CH0314K2H	27 mm	31 mm	32 mm	10	416.00	482.00
CH0314L1H	CH0314L2H	29 mm	33 mm	32 mm	10	458.00	519.00
CH0314M1H	CH0314M2H	31 mm	36 mm	35 mm	10	660.00	733.00
CH0314N1H	CH0314N2H	33 mm	38 mm	38 mm	10	807.00	880.00
CH0314O1H	CH0314O2H	35 mm	45 mm	36 mm	10	959.00	996.00
CH0314P1H	CH0314P2H	38 mm	42 mm	40 mm	10	1741.00	1179.00
CH0314Q1H	CH0314Q2H	40 mm	49 mm	40 mm	10	1106.00	1540.00
CH0314R1H	CH0314R2H	44 mm	58 mm	45 mm	10	1418.00	1863.00
CH0314S1H	CH0314S2H	48 mm	54 mm	33 mm	10	2762.00	2884.00
CH0314T1H	CH0314T2H	51 mm	60 mm	57 mm	10	3000.00	3122.00
CH0314X1H	CH0314X2H	75 mm	90 mm	39 mm	10	5109.00	5231.00





Rubber Stoppers - One Hole, Pack of 50

- Made from superior rubber, 1 hole
- Resistant to acids, alkalis and ammonia
- 7 pieces each 8, 9, 11, 13, 15 mm bottom dia.
- 3 pieces each 17, 21, 23, 25, 29 mm bottom dia.
- Pack of 50

Code	Specs	Pack Size	Price/Pack ₹
CH0316 ①	Assorted Rubber Stoppers - One Hole	50	1120.00



Rubber Stopper - Two Hole, Pack of 50

- Made from superior rubber, 2 hole
- Resistant to acids, alkalis and ammonia
- 10 pieces each 13, 15, 17, 18 mm bottom dia.
- 3 pieces each 21, 23 mm bottom dia.
- 2 pieces each 25, 29 mm bottom dia.
- Pack of 50

Code	Specs	Pack Size	Price/Pack ₹
CH0317 ②	Assorted Rubber Stoppers - Two Hole	50	1250.00



Silicone Stoppers - Solid ③

- Made from pure silicon rubber
- Safe for use upto 430°F / 225°C
- Available in following sizes

Code	Bottom	Top	Length	Pack Size	Price/Pack ₹
CH0319A	6 mm	8 mm	16 mm	20	373.00
CH0319B	8 mm	10.5 mm	20 mm	20	384.00
CH0319C	9 mm	11.5 mm	20 mm	20	396.00
CH0319D	10 mm	12.5 mm	20 mm	20	528.00
CH0319E	11 mm	14 mm	24 mm	20	616.00
CH0319F	13 mm	16 mm	24 mm	20	791.00
CH0319G	15 mm	18 mm	24 mm	20	1011.00
CH0319H	17 mm	20 mm	26 mm	20	1231.00
CH0319I	18 mm	21 mm	26 mm	20	1408.00
CH0319J	19 mm	22 mm	28 mm	20	1628.00
CH0319K	21 mm	24 mm	28 mm	20	1980.00
CH0319L	23 mm	26 mm	28 mm	20	2376.00
CH0319M	25 mm	28 mm	28 mm	10	2816.00
CH0319N	27 mm	31 mm	32 mm	10	3300.00
CH0319O	29 mm	31 mm	32 mm	10	4268.00
CH0319P	31 mm	36 mm	35 mm	10	5500.00
CH0319Q	33 mm	38 mm	38 mm	10	7700.00
CH0319R	35 mm	45 mm	36 mm	10	6967.00
CH0319S	38 mm	42 mm	40 mm	10	9167.00
CH0319T	40 mm	49 mm	40 mm	10	11148.00



EISCO TUBING

- Research grade for industrial and academic research labs
- Resistant to acids and alkalis
- Neoprene, nitrile, rubber, extra soft and high pressure varieties available

Rubber Tubing 1

- Soft quality, resistant to acid and alkalis
- Shore hardness: 45A
- Tolerance: $\pm 0.2\text{mm}$
- Temp.: -20°C to $+70^{\circ}\text{C}$

Code	Bore	Wall Thickness	Pack Size	Price/Pack ₹
CH0761A01	3 mm	1.0 mm	10 Mtr	508.00
CH0761A	4 mm	1.0 mm	10 Mtr	591.00
CH0761B	5 mm	1.5 mm	10 Mtr	738.00
CH0761C	6 mm	1.5 mm	10 Mtr	956.00
CH0761C01	6.5 mm	1.5 mm	10 Mtr	1179.00
CH0761D	7 mm	1.5 mm	10 Mtr	1032.00
CH0761E	8 mm	2.0 mm	10 Mtr	1357.00
CH0761F	9 mm	2.0 mm	10 Mtr	1472.00
CH0761G	10 mm	2.0 mm	10 Mtr	1802.00
CH0761H	12 mm	2.5 mm	10 Mtr	2273.00
CH0761I	15 mm	3.0 mm	10 Mtr	3862.00
CH0761K	15 mm	3.0 mm	10 Mtr	3911.00
CH0761L	16 mm	3.0 mm	10 Mtr	4094.00
CH0761M	19 mm	3.0 mm	10 Mtr	4620.00
CH0761N	20 mm	3.5 mm	10 Mtr	4718.00
CH0761O	25 mm	3.5 mm	10 Mtr	5622.00



1



Rubber Cones

- Pack of 7 rubber cones for use with buchner funnels of various sizes

Code	Capacity	Pack Size	Price/Pack ₹
CH0421 ¹	Rubber cones	7	649.00



Rubber Tubing - High Pressure

- High pressure, extra soft quality rubber tubing
- Shore Hardness: 50-55A
- Tolerance: ±0.2mm
- Temp.: -20°C to +70°C

Code	Bore	Wall Thickness	Pack Size	Price/Pack ₹
CH0762A	3 mm	4.5 mm	10 Mtr	1467.00
CH0762B	4.5 mm	4.5 mm	10 Mtr	1680.00
CH0762C	6 mm	4.5 mm	10 Mtr	1964.00
CH0762D	8 mm	4.5 mm	10 Mtr	2316.00
CH0762E	9 mm	4.5 mm	10 Mtr	2603.00
CH0762E01	10 mm	6.0 mm	10 Mtr	4119.00
CH0762F	12.5 mm	6.5 mm	10 Mtr	4308.00



Make sure to sign up for our **monthly** email



Silicone Tubing ¹

- Seamless, non-toxic, translucent and very flexible
- Typically used for connecting laboratory glassware
- Suitable for use over a temperature range of -50°C to 250°C
- Meets USP VI standards and FDA 21 CFR 177.2600 compliant
- Shore Hardness: 45A

Code	Bore	Wall Thickness	Pack Size	Price/Pack ₹
CH0764A	2 mm	1.0 mm	10 Mtr	1051.00
CH0764B	3 mm	1.0 mm	10 Mtr	1259.00
CH0764B01	3 mm	1.5 mm	10 Mtr	1802.00
CH0764D	4 mm	1.0 mm	10 Mtr	1784.00
CH0764E	4 mm	1.5 mm	10 Mtr	1980.00
CH0764F	5 mm	1.0 mm	10 Mtr	1980.00
CH0764G	5 mm	1.5 mm	10 Mtr	2616.00
CH0764I01	6 mm	1.0 mm	10 Mtr	3324.00
CH0764I	6 mm	1.5 mm	10 Mtr	2616.00
CH0764J	6 mm	2.0 mm	10 Mtr	4033.00
CH0764J01	6 mm	3.0 mm	10 Mtr	6500.00
CH0764M01	7 mm	1.5 mm	10 Mtr	3333.00
CH0764M	7 mm	2.0 mm	10 Mtr	4222.00
CH0764O01	8 mm	1.5 mm	10 Mtr	3678.00
CH0764O	8 mm	2.0 mm	10 Mtr	4200.00
CH0764P	8 mm	3.0 mm	10 Mtr	6056.00
CH0764Q	9 mm	2.0 mm	10 Mtr	5000.00
CH0764R	10 mm	2.0 mm	10 Mtr	6056.00
CH0764S	10 mm	3.0 mm	10 Mtr	9067.00
CH0764T01	12 mm	2.0 mm	10 Mtr	6400.00
CH0764T	12 mm	3.5 mm	10 Mtr	10978.00
CH0764V	14 mm	3.0 mm	10 Mtr	11833.00
CH0764W	15 mm	4.0 mm	10 Mtr	17333.00
CH0764X	16 mm	3.0 mm	10 Mtr	14458.00
CH0764X01	25 mm	4.0 mm	10 Mtr	19911.00

